

REFERENCE GUIDE

Adapters and Indices for the G4 Sequencing Platform

Singular adapters and unique dual indices maximize flexibility for the independently addressable lanes of the G4 Sequencing Platform.

The G4™ Sequencing Platform by Singular Genomics™ is the most powerful benchtop sequencer to date. The G4 Sequencing Platform and Singular sequencing chemistry have been designed from the ground up, resulting in increased speed and flexibility, and support a wide range of library preparation workflows. Singular Genomics has developed an adapter strategy with 96 unique dual indices that enable efficient use of each independently addressable lane and flow cell, and allows for flexibility in the index sequence lengths. This reference guide describes the adapters and indices, explains how to integrate them in library preparation workflows, and provides guidelines for using custom libraries on the G4 Sequencing Platform.

SINGULAR ADAPTERS FOR LIBRARY PREP

Singular Adapter Sequences

Singular adapters are used during library preparation to add several functional nucleotide tags to each end of the inserts to be sequenced (Fig 1). At the 5' ends, Singular proprietary platform sequences S1 and S2 are attached as anchors for the formation of clusters on the flow cell. The SP1 and SP2 tags, which are identical to the SP1 and SP2 sequencing primers used in many existing applications, are positioned directly adjacent to the insert. When constructing libraries for multiplexed reads, index 1 and index 2 sequences are placed in between S1 and SP1, and S2 and SP2 respectively. The inserts and indices are sequenced using the primers described in Table 1.




Figure 1: Singular Library with Dual Indices. **A.** Structure of Singular library, arrows denote direction of sequencing. **B.** Final sequence of adapters on a library with dual indices.

Table 1: Sequencing Primers.

Primer	Sequence Order	Sequenced Segment	Sequence of Primer
S1	1	Index 1	5'-ACAAAGGCAGCCACGCACTCCTTCCCTGT
SP1	2	Insert Read 1	5'-AACTCTTTCCCTACACGACGCTCTTCCGATCT
S2	3	Index 2	5'-CTCCAGCGAGATGACCTCACCAACCACT
SP2	4	Insert Read 2	5'-GTGACTGGAGTTCAGACGTGTGCTCTTCCGATCT

Universal Stem-Loop Adapter

Singular Genomics uses a universal stem-loop adapter during library preparation as first step to introduce the nucleotide sequences necessary for cluster generation and sequencing. The stem-loop adapter consists of short complementary sequences at the 5' and 3' ends, a non-complementary segment in the center with uracil in the middle, and a T-overhang at the 3' end (Figure 2). The 3' end contains the SP1 sequencing



Figure 2: Sequence of the Singular Universal Stem-Loop Adapter.

to the 3' ends (Fig. 3, step 1). This A-tail provides an efficient target to ligate with the 3' T-tail of the stem-loop adapter (Fig. 3, step 2). Next, the cleavable site in the middle of the loop is excised (Fig. 3, step 3) using a Uracil-DNA Glycosylase¹ (UDG) in appropriate buffer conditions or a suitable alternative. Cleaving the backbone opens the loops and results in a linear fragment with the SP1 sequence at the 5' end, and the SP2' sequence at the 3' end. These sequences are then used to add the Singular S1 cluster sequence and index 1 at the SP1 end, and Singular S2 cluster sequence and index 2 at the SP2 end through a limited amplification step (Fig. 3, step 4-6). Primers without indices are used if no multiplexing is required (Table 2).

primer sequence, while the 5' end is complementary to the standard SP2 sequencing primer (SP2'), without the A at the 5' end. SP1 and SP2 are sequencing primers widely used in existing sequencing protocols.

Standard Library Preparation

Adding the necessary sequences for cluster generation and sequencing to your library is a straight-forward procedure.

After you process and fragment a DNA sample, add an A-tail




Figure 3: Singular Standard Library Preparation.

Primer	Sequence of Primer
S1-Index-SP1	5'- ^{S1} ACAAAGGCAGCCACGCACTCCTTCCTGT[index 1] ^{SP1} ACACTCTTTCCTACACGACGCTTCCGATCT
S2-Index-SP2	5'- ^{S2} CTCCAGCGAGATGACCCTCACCAACCACT[index 2] ^{SP2} GTGACTGGAGTTCAGACGTGTGCTCTTCCGATCT
S1-SP1	5'- ^{S1} ACAAAGGCAGCCACGCACTCCTTCCTGTACACTCTTTCCTACACGACGCTTCCGATCT
S2-SP2	5'- ^{S2} CTCCAGCGAGATGACCCTCACCAACCACTGTGACTGGAGTTCAGACGTGTGCTCTTCCGATCT

Table 2: PCR primers for introducing cluster sequences (S1 and S2) with or without indices.

The steps described above result in a library that is ready for clustering and sequencing on the G4 Sequencing Platform (Fig.1). The S1 and S2 sequences are used to form clusters on the flow cell. Then, depending on the application, the inserts and indices are sequenced using the primers described in Table 1.

PCR-Free Library Preparation

PCR-free library preparation has certain advantages, including eliminating PCR bias and a simplified workflow, though with lower yield. For PCR-free sequencing, Singular uses a different stem-loop adapter. In addition to the uracil and SP1 and SP2' sequences, the PCR-free stem-loop adapter also contains the S1 and the complementary sequence to S2 (Fig.4A). For indexed applications, Singular has embedded dual indices within the adapter design (Fig.4B). This means no amplification is required to add the index sequences, which simplifies the workflow (Fig. 4C).




Figure 4: Singular PCR-free Library Preparation. **A.** Sequence of Singular stem-loop adapter for PCR-free library preparation. **B.** Like A, but with dual indices. **C.** Singular PCR-free library preparation steps.

SINGULAR INDEXING STRATEGY

Indices are short DNA barcodes added during library prep that allow you to correctly assign reads to their corresponding sample via demultiplexing. When pooling multiple libraries within a single lane, you should always assign a unique combination of indices to each library. An ideal index design avoids misassociation of reads, and the assignment process is resistant to errors, including those introduced during index synthesis. The index sequences should be as divergent as possible while having similar GC content and avoiding extended homopolymer stretches. Indices that are many changes (insertions, deletions, and/or substitutions) away from other index sequences minimize the possibility of misassociation.

Singular dual indices have been designed with these principles in mind, and have the following features:

- Singular Genomics provides 96 unique dual indices (UDIs) of 12 bases for sample indexing.
- Robust UDI design allows for efficient demultiplexing with as few as 8 cycles, which enables efficient demultiplexing of 96 libraries. Extending reads to 12 cycles assures maximal differentiation between samples.
- The indices have been empirically tested to ensure even library representation in multiplex experiments.
- Quality control analysis of indices ensures high purity oligomer synthesis.
- Rapid demultiplexing can be performed on or off instrument via the Singular Genomics Demultiplex Tool, an accelerated version of fgbio's DemuxFastqs.

The complete set of Singular UDIs is listed in Appendix 1: Singular Dual Indices.

It is best practice to use dual indices, which enable identification and removal of mis-associated indices that can arise during library prep and index hopping. While index hopping in cluster amplification is rare (< 0.1% of reads), levels close to 2% and up to 10% have been reported for other platforms^{2,3,4}. We have measured index hopping on the G4 Sequencing Platform, and have observed very low levels of index hopping in the order of 0.07% (for details see Appendix 2: Index Hopping), indicating index hopping is not a substantial source of error for the G4 Sequencing Platform.

ADAPT LIBRARY FORMAT TO SINGULAR SEQUENCING

Considerations When Using Custom Sequencing Primers

If you design your own custom adapters and primers, the adapters need to contain S1 and S2 sequences for clustering on the G4. The S1 and S2 sequences need to be located on the ends of the fragment, so your custom sequences and insert will be amplified during cluster generation (Fig. 5). Note that you must validate custom primers yourself; Singular Genomics cannot guarantee performance or compatibility of custom primers.




Figure 5: Outline of Custom Library during Cluster Generation. Indicated are the custom parts of the adapter, as well as the Singular proprietary sequences needed for cluster generation with Singular S1 and S2 oligos.

Your custom sequencing primers need to account for the melting temperature (T_m) relative to sequencing temperature. Also be aware of the salt concentration of the buffer for the sequencing primers, as it impacts hybridization efficiency and may impact amplification and sequencing efficiency. Custom sequencing primers need to be able to bind at 66°C (while assuming 50 mM Na^+). Ideally, your custom primers are very similar in GC content and length to SP1 and SP2 sequencing primers (52% and 33 bp). Also avoid short GC-rich sequencing primers. Contact the Singular Genomics Customer Care team for additional technical guidance.

Existing Library Conversion

Converting existing libraries to Singular libraries is usually straightforward. Many libraries already contain SP1 and SP2 sequences adjacent to the insert. You can do a simple PCR using either the S1-index1-SP1 and S2-index2-SP2 primers, or, if no index is required, use the S1-SP1 and S2-SP2 primers (Table 2 and Fig. 6). After this, use a PCR purification kit to clean up your product and follow the regular workflow for clustering on the G4.




Figure 6: Strategy for Adapting Existing Libraries to Singular Sequencing.

Note that the SP1 and SP2 sequencing primer sequences the G4 uses are identical to the ones used in many legacy platforms, while the indices and S1 and S2 are unique for Singular Genomics.

If you have an existing library without SP1 and SP2 sequencing primers, you will have to design your own conversion primers. They should minimally consist of the S1 and S2 sequences at the 5' end for clustering on the G4, and your custom sequencing primer sequences for anchoring on your library fragment (see Design Custom Libraries for additional guidelines, and contact Singular Genomics Customer Care for help).

Adapter Trimming

Trimming adapters from Singular Genomics data can be done using common NGS bioinformatic tools such as cutadapt⁵ or BBDuk⁶. We recommend using at least 12 bp of the adapter sequence to trim. For standard runs using the SP1 and SP2 sequencing primers, the sequences that can be used to trim R1 and R2 are listed in Table 3.

Table 3: Trimming Sequences.

Sequencing Primer	Sequenced Segment	Recommended Trimming Sequence
SP1	Read 1	AGATCGGAAGAGCACA
SP2	Read 2	AGATCGGAAGAGCGTC

RESOURCES

Customer Care

For any questions regarding this document, or for further assistance, contact Singular Genomics Customer Care:

Email: care@singulargenomics.com

Call: +1-442-SG-CARES (1-442-742-2737)

Website: www.singulargenomics.com

Ordering Products

To order the products described in this reference guide from Singular Genomics, use the following part numbers:

Part Number	Product
700110	SG UDI - Set of 96
700111	SG UDI - Set of 24 [Set A]
700112	SG UDI - Set of 24 [Set B]
700115	G4 Non-Indexed PCR-Free Library Prep Adapters (24 Reactions)
700116	G4 Indexed PCR-Free Library Prep Adapters (24 Reactions)
700117	G4 Universal Library Prep Adapters (96 Reactions)
700118	G4 Universal Library Prep Adapters (24 Reactions)
700119	G4 Non-Indexed Library Prep Primers (24 reactions)

REFERENCES

1. Pearl, LH. Structure and function in the uracil-DNA glycosylase superfamily. *Mutat Res.* 2000; 460(3-4):165-81. doi: 10.1016/s0921-8777(00)00025-2
2. Costello M, Fleharty M, Abreu J, Farjoun Y, Ferriera S, Holmes L, Granger B, Green L, Howd T, Mason T, Vicente G, Dasilva M, Brodeur W, DeSmet T, Dodge S, Lennon NJ, Gabriel S. Characterization and remediation of sample index swaps by non-redundant dual indexing on massively parallel sequencing platforms. *BMC Genomics.* 2018; 19:332. doi: 10.1186/s12864-018-4703-0
3. Effects of Index Misassignment on Multiplexing and Downstream Analysis. index-hopping-white-paper-770-2017-004.pdf
4. Sinha R, Stanley G, Gulati GS, Ezran C, Travaglini KJ, Wei E, Chan CKF, Nabhan AN, Su T, Morganti RM, Conley SD, Chaib H, Red-Horse K, Longaker MT, Snyder MP, Krasnow MA, Weissman IL. Index switching causes "spreading-of-signal" among multiplexed samples in Illumina HiSeq 4000 DNA sequencing. *bioRxiv.* 2017; 125724. doi: 10.1101/125724
5. <https://cutadapt.readthedocs.io/en/stable/>
6. <https://jgi.doe.gov/data-and-tools/software-tools/bbtools/bb-tools-user-guide/bbduk-guide/>
7. Schwendinger J. Sequence with confidence: understand index hopping and how to resolve it. <https://www.10xgenomics.com/blog/sequence-with-confidence-understand-index-hopping-and-how-to-resolve-it>
8. Farouni R, Djambazian H, Ferri LE, Ragoussis J, Najafabadi HS. Model-based analysis of sample index hopping reveals its widespread artifacts in multiplexed single-cell RNA-sequencing. *Nature Communications.* 2020; 11: 2704
9. Van der Valk T, Vezzi F, Ormestad M, Dalén L, Guschanski K. Index hopping on the Illumina HiSeqX platform and its consequences for ancient DNA studies. *Mol Ecol Resour.* 2020; 20(5):1171-1181. doi: 10.1111/1755-0998.13009
10. Ros-Freixedes R, Battagin M, Johnsson M, Gorjanc G, Mileham AJ, Rounsley SD, Hickey JM. Impact of index hopping and bias towards the reference allele on accuracy of genotype calls from low-coverage sequencing. *Genet Sel Evol.* 2018; 50(1):64. doi: 10.1186/s12711-018-0436-4.

APPENDIX 1: SINGULAR DUAL INDEX SEQUENCES

Indicated are the Singular dual index sequences in each set. S1 and S2 in the index name indicate the index can be sequenced using the S1 primer and S2 primer, respectively. After each 12-base index sequence, there is a single nucleotide spacer before SP1 or SP2 sequences start. See also Table 2.

Set A

Dual Index 1	Dual Index 2
S1-Index1 TAAGACCCTACT T	S2-Index1 GGGACATATTGA C
S1-Index2 CGAAGTACATCC C	S2-Index2 TAGGACGTAACG G
S1-Index3 TAGCCTTCCAAA C	S2-Index3 AGTATGGCAAGA C
S1-Index4 GCCTTTCAAGTC C	S2-Index4 TAGAGTCGTCGT G
S1-Index5 CAACGGTTCCGG C	S2-Index5 ACGTTTCGCTCG T
S1-Index6 GTTGCATGGCCC T	S2-Index6 TAGGGAACGATG G
S1-Index7 ATCGTTGCTATC G	S2-Index7 ATGACTCCGCAT T
S1-Index8 CCTCGAATTTCAT C	S2-Index8 GGTTGCTACCGG T
S1-Index9 TGAACGTCCGCC C	S2-Index9 TCCTCGATTGAA G
S1-Index10 CATCTAGCAAGC T	S2-Index10 ATGTAGCGTCTC A
S1-Index11 TATCGAGGCAAC C	S2-Index11 CATCATGCGTAC C
S1-Index12 GAGACGTAGCAA G	S2-Index12 ACCTTGACCGGG G
S1-Index13 ATCATGCGCCCG T	S2-Index13 TTGACGAGATCT A
S1-Index14 AGGAGCTAGGGA G	S2-Index14 GGGCTAATGTCA T
S1-Index15 ATCGACCATGCT C	S2-Index15 TTAGGAGCGAAC T
S1-Index16 TGCGAATCGACA C	S2-Index16 GTACATCGAGTA G
S1-Index17 ATGTTCCCCTCT T	S2-Index17 AGGCTTTGTCAT T
S1-Index18 TCGCTCATCTAG A	S2-Index18 CGACGATATTTG G
S1-Index19 CCTAAGGTAAAC C	S2-Index19 AAACCTCCGTTGT T
S1-Index20 GAATAGCGCTTA T	S2-Index20 ACGTACCAAGAC C
S1-Index21 CGATGTACATCC T	S2-Index21 TCCGATGTCGGC G
S1-Index22 CAAGTCGAAACC T	S2-Index22 AGGTTACCGCGT A
S1-Index23 GTAACGGATAGC T	S2-Index23 ATGCCGAAACGT T
S1-Index24 GAAGCTTGGTCA G	S2-Index24 ACATACGCGGGG G

Set C

Dual Index 1	Dual Index 2
S1-Index49 TCAAGAGCGGAG G	S2-Index49 CGGACTTTTGTA G
S1-Index50 TTACCCGTAGAA A	S2-Index50 TCCATTGCTTCT A
S1-Index51 GCTCTCAATCGG G	S2-Index51 GTCTCATGGCGG A
S1-Index52 GTCTACGTTTAC C	S2-Index52 TGACTTGAGAA A
S1-Index53 TCCGTATGAGAC C	S2-Index53 ACCGTATCCGAT C
S1-Index54 CGCCAATACGTC A	S2-Index54 CTATGGGACGGT T
S1-Index55 GATGGTCTAGCA C	S2-Index55 TAGTTCCCATTC C
S1-Index56 CTCGCTTAAGGC G	S2-Index56 CTCCAAGACATC A
S1-Index57 GGCAACATGGGT C	S2-Index57 ATTACCGCGGTA C
S1-Index58 AGACTCTCATCA A	S2-Index58 CTAAGCTCCTAA C
S1-Index59 TGACAAGGTCAA A	S2-Index59 TGAACGTCCCTC T
S1-Index60 CGGTATGTCATC G	S2-Index60 CCATGCAATCTA C
S1-Index61 GACTCATGAATG A	S2-Index61 TAGCGTTTCATTG G
S1-Index62 CGTAGACATTGA T	S2-Index62 AGCATTCGAGCC T
S1-Index63 CATTGCTCCCT G	S2-Index63 AACCTCGAACAT A
S1-Index64 ACAATCGGGGAC T	S2-Index64 CAAATCGCGGAA A
S1-Index65 GGACTTAGAGCG C	S2-Index65 GTCAAGAGGTTA A
S1-Index66 GACCGATTCTCG A	S2-Index66 TTAAGGCCGGGA G
S1-Index67 TGGAACCCGAG T	S2-Index67 TCGAAAGGGAAA C
S1-Index68 GGCCTAATGGAA G	S2-Index68 GGAGTCAAATAG T
S1-Index69 TTGTACGCGTAC A	S2-Index69 CCGTTCTATACA A
S1-Index70 ATGTCGAGTTGC A	S2-Index70 TGCGAATGCGAA C
S1-Index71 CTTCTGACCTCC A	S2-Index71 GCGATCATGACT C
S1-Index72 TTAGGTCGAGAA A	S2-Index72 CCGAAATCCAAC C

Set B

Dual Index 1	Dual Index 2
S1-Index25 AACCCGTAACCA T	S2-Index25 GGATCTAGGACG G
S1-Index26 AATGCTCCCTA C	S2-Index26 TCGACTCTCCGT T
S1-Index27 GTATGACGGATG G	S2-Index27 TACGCTAGACAA A
S1-Index28 GCAAAGCTTGGA C	S2-Index28 ATTTGCGGTAAG T
S1-Index29 TCTAACCGGCTA G	S2-Index29 CTAGCCAACGCC G
S1-Index30 ATTGGAGCCCGC C	S2-Index30 TGATAGCCGGTT C
S1-Index31 TGTCCGATCTAT A	S2-Index31 GCATGGTTCCTA T
S1-Index32 ACATCGCATGTT C	S2-Index32 GACGGTAATGAG T
S1-Index33 ACTTCCGACAAT A	S2-Index33 CGAAGTCATCAA C
S1-Index34 GCTCCTATGCCT T	S2-Index34 TGGACTACAAAC T
S1-Index35 CATAACGCGAAT G	S2-Index35 ATACGCGACATT G
S1-Index36 CGGACAATGATT G	S2-Index36 TCCCTATGATAC C
S1-Index37 GCTACTTGGA A	S2-Index37 GAGCTCTACGCA T
S1-Index38 TAAGCGAGTAGT C	S2-Index38 CTAGAGGTTACC T
S1-Index39 TGCTTGAATCCG G	S2-Index39 AACCCCTGTTTCG T
S1-Index40 GACATCGTCGGG C	S2-Index40 GTCGTAAGGGGT T
S1-Index41 AGCTATGGGACG A	S2-Index41 TATACCCGGCCC A
S1-Index42 ACGACTAGGCTC T	S2-Index42 TGCTAAGCGAGC C
S1-Index43 TACCGTACGATA A	S2-Index43 CGTCCTAAAAC T
S1-Index44 TAGAAGGCGCGT C	S2-Index44 ATATCGGGTAAG A
S1-Index45 CGTCATCAAGGA G	S2-Index45 GAGACGTCTCTA A
S1-Index46 CCGTAAGATAGA C	S2-Index46 GCTTAACGATCA T
S1-Index47 TAGACTCGTTTC C	S2-Index47 CGCTAGTACTAT G
S1-Index48 TATCGGCTTGGT C	S2-Index48 ACGGGTTATTAG A

Set D

Dual Index 1	Dual Index 2
S1-Index73 CTAGCTCTTCGT T	S2-Index73 TCGGAGTTTAC C
S1-Index74 CTTGTCCAACCT G	S2-Index74 CTATCCGTCCCG C
S1-Index75 CTTAGCGACCCA T	S2-Index75 TAACGCGTACCC C
S1-Index76 CGTAGGTAAACA A	S2-Index76 CTCGTACTTAGC A
S1-Index77 AGCATTCATGT G	S2-Index77 GACTCGAAATGA A
S1-Index78 TCGTTACCAACG C	S2-Index78 TTAACCGGCCGA C
S1-Index79 TTGCTAGGACAT G	S2-Index79 TTCTAGCAACC C
S1-Index80 CGAGACTTCTAC A	S2-Index80 GAATAGCGTCCC T
S1-Index81 GGTCTATGTTTG T	S2-Index81 AATGACGGATGT G
S1-Index82 GATGCCATAGTA G	S2-Index82 AGTAGCTCGTCC C
S1-Index83 GTACGAGTTCCCT C	S2-Index83 TGCCATTCTCC T
S1-Index84 TTCCATCGGTAG A	S2-Index84 ACGCCATAACCG T
S1-Index85 ACGTATCATCT T	S2-Index85 CGTTAACCCGCT T
S1-Index86 GTCCAAGAGTTC C	S2-Index86 TTGAGCCTTGCT A
S1-Index87 CGAATGGTAAGA T	S2-Index87 ACCCGTAACAAG G
S1-Index88 AGGTTTCGGTAT C	S2-Index88 AAAGGTCGCCCC A
S1-Index89 GGTTCAAGATAC T	S2-Index89 TACGAGCTCCTC A
S1-Index90 GTTAAGCGGGCG A	S2-Index90 GGCATAGAGCGA T
S1-Index91 AAGCTACGCGAA T	S2-Index91 GGTAGCATTCGG T
S1-Index92 TTGAGCGTCCGA C	S2-Index92 CGACAAGTTCGA G
S1-Index93 CCTATGCAATTG T	S2-Index93 GGCAATGTACTT C
S1-Index94 CTATTGCCTACG A	S2-Index94 CCGGTAAACCCA C
S1-Index95 CCTCTTGAAGAG A	S2-Index95 GCCCATTAGATT G
S1-Index96 TGGGTTACGGGC G	S2-Index96 CAAGGTTCTCTT C

APPENDIX 2: INDEX HOPPING

Indices are widely used for multiplexing to identify which reads belong to which samples, but sometimes errors happen. There are several causes for index misassociation. Some of them, like poor oligo synthesis or well-to-well contamination, can be eliminated through good practices. Others are caused by the chemistry employed during cluster formation (for example, index hopping).

A common cause of index hopping is contamination with free adapters and primers during cluster generation. These free adapters can hybridize with the complementary sequence of an adapter from a different sample and extend. This free hybrid molecule can seed a cluster with an incorrect index for that insert in a different well. There are other causes of index hopping, like low diversity in sequence that allows two fragments to bind to each other. These fragments can then form hybrids that extend and seed clusters with mismatched indices.^{2,7}

While index hopping in bridge amplification usually stays well below 1% of the reads, levels close to 2%, with levels as high as 10% have been reported on other SBS platforms.^{2,3,4} These levels can lead to serious problems in downstream analysis in a variety of applications, like somatic variant detection and differential expression.^{8,9,10} We therefore wanted to characterize the potential for index hopping on the patterned flow cells of the G4 Sequencing Platform.

We ran 4 lanes each with 24 indexed *Salmonella* whole-genome samples using different unique dual indices for each lane. The data were analyzed to detect any unexpected dual index pairs between the 24 dual indices, which are indicative of index hopping events. As shown in Fig. 7, in all lanes, index hopping events were rare, with on average only 0.07% of reads having aberrant dual index combinations.




Figure 7: Counts of dual index pairs between 24 dual indices on each lane. Expected combinations are located on the diagonals.

Even though index hopping on the G4 Sequencing Platform occurs at low frequency, its impact on downstream analysis can still be minimized further. Singular Genomics uses unique dual indices and filters out unexpected dual index pairs during demultiplexing. In addition, there are laboratory best practices to prevent index misassignment:

- Make sure to completely remove free adapters after the ligation step.
- Store prepared libraries at the recommended temperature of -20° C.
- Employ good practices to prevent cross-contamination of oligos during library preparation.
- Singular Genomics provides adapters of high purity. If you plan to use custom adapters, make sure they are of the highest quality.

In conclusion, a low level of inherent index hopping, dual indices, and best practices combine to ensure that index hopping is not a substantial source of artifacts for the G4 Sequencing Platform.



SINGULAR
GENOMICS

For Research Use Only.
Not for use in diagnostic procedures.

CUSTOMER CARE

Our Singular focus is you.

Email: care@singulargenomics.com

Call: +1-442-SG-CARES (1-442-742-2737)

Website: www.singulargenomics.com



PN 600007 Rev. 0, June 29, 2022